



**NANYANG  
TECHNOLOGICAL  
UNIVERSITY**  
SINGAPORE

# MS4672 – Week 5

## High Throughput Experimental Methods for Materials Discovery

Emha Bayu Miftahullatif (EMHA0001@e.ntu.edu.sg)  
Pethe Shreyas Dinesh (SHREYASD002@e.ntu.edu.sg)

Also thanks to:  
Dr. Danny Ren Zekun

School of Materials Science and Engineering



# Today's class

- 1st half : Introduction to High-Throughput Experimentation (120 mins)
  - Introduction to High Throughput experimentation (20 mins)
  - Opentrons automated pipetting robots + Code walkthrough (40 mins)
  - In-class project: automated solution mixing (60 mins)
- 2<sup>nd</sup> half: OpenCV image manipulation (60 mins)
  - Introduction to OpenCV (30 mins)
  - Code walkthrough (20 mins)
  - Homework (10 mins)

# **1st half**

## **Introduction to High-Throughput Experimentation**



# Conventional Discovery

## ① Serendipity-driven

Gorilla glass (1952)



*Pure luck*

## ② Persistent trial and error


Light bulb filament (1880)




*Hypothesis-driven*  
*Slow and expensive*

# Why?

Because materials synthesis is complex




Even more complex to make device out of it



The new mainstream:

# High Throughput Experimentation + self-driving lab




- Accelerated data acquisition
- Enhanced Reproducibility
- Resource and Cost Efficiency
- AI-Driven Optimization
- Continuous Operation
- Broad Parameter Mapping
- Seamless Data Integration



The new mainstream:

# High Throughput Experimentation + self-driving lab








# Types of modules to accomplish the tasks


## Modules for making things out of stuff




## Modules for Testing



## Modules for moving stuff



## Decision-making brain



# Types of experimental techniques in HTE


- Point-wise experiments (e.g. mechanical, spectrum analysis, electrical etc...)
- Parallel/Simultaneous experiments (e.g. impedance spectroscopy, gas adsorption etc...)
- Proxy experiments (Optical, electrical measurements)
- Ultrafast analysis



Nanoindenter




Potentiostat




Hyperspectrum

High-resolution scans, 100x faster than conventional CT



18650 BATTERY CELL  
10 second Ultra-Fast CT scan



PLASTIC BOTTLE  
24 second Ultra-Fast CT scan

Ultrafast CT

<https://www.bruker.com/en/products-and-solutions/test-and-measurement/nanomechanical-test-systems/hysitron-ti-980-nanoindenter.html>

[https://www.biologic.net/product\\_category/potentiostats-galvanostats/](https://www.biologic.net/product_category/potentiostats-galvanostats/)


<https://www.lumafield.com/product-feature/ultra-fast-ct-scanning>

Wikipedia




# Workflow example

## Optimizing perovskite precursors for printable solar cells





Next:  
We will learn to use this



Miftahullatif, Emha B.; et al. Machine-Learning-Driven in-Device Optimization of All-Printed Perovskite Solar Cells. *ACS Energy Lett.* **2025**, 10 (8), 3952–3961.

# Introduction to Opentrons

- 3-axis Liquid handling robot
- Picks liquid from reservoir (resource), pipettes into receptacle, disposes tips in trash
- Python controlled
- Standardized labware



# Introduction to Opentrons



## Modules



Pipettes



thermocycles



Heater



Heater-shaker

## Labwares




## Protocol

```
1 protocol_modified.py X
2
3 from opentrons import protocol_api
4
5 # Protocol metadata - required for all protocols
6 metadata = {
7     'protocolName': 'Basic Liquid Handling',
8     'author': 'Student Name',
9     'description': 'Simple protocol for learning OT-2 operation',
10    'apiLevel': '2.13'
11}
12
13 def run(protocol: protocol_api.ProtocolContext):
14     """
15     This function defines what the robot will do.
16     It's called automatically when the protocol runs.
17     """
18
19     # 1. Load labware (define what equipment is on the deck)
20     # Slot numbers 1-11 refer to positions on the robot deck
21     plate = protocol.load_labware('corning_96_wellplate_384ul_flat', '1')
22     tiprack = protocol.load_labware('opentrons_96_tiprack_300ul', '11')
23     reservoir = protocol.load_labware('opentrons_10_tuberack_falcon_4x50ml_6x15ml_conical', '3')
24
25     # 2. Load pipette (define which pipette to use)
26     # 'right' or 'left' refers to the pipette mount position
27     pipette = protocol.load_instrument(
28         'p300_single_gen2',
29         'left',
30         tip_racks=[tiprack]
31     )
32
33     # 3. Protocol steps (define the liquid handling operations)
34
35     # Pick up a tip
36     # pipette.pick_up_tip()
37     pipette.pick_up_tip(tiprack['A1'])
38
39     # Aspirate 20ul from reservoir well A1
40     pipette.aspirate(20, reservoir['A3'])
41
42     # Dispense 20ul to plate well A1
43     pipette.dispense(20, plate['A1'])
44
45     # Repeat for multiple wells
46     for well in ['A2', 'A3', 'A4', 'A5', 'A6']:
47         pipette.aspirate(20, reservoir['A3'])
48         pipette.dispense(20, plate[well])
49
50     # Drop the tip when finished
51     pipette.drop_tip()
```



# Protocol

```
protocol_modified.py X
protocol_modified.py > ...
1  from opentrons import protocol_api
2
3  # Protocol metadata - required for all protocols
4  metadata = {
5      'protocolName': 'Basic Liquid Handling',
6      'author': 'Student Name',
7      'description': 'Simple protocol for learning OT-2 operation',
8      'apiLevel': '2.13'
9  }
10
11  def run(protocol: protocol_api.ProtocolContext):
12      """
13      This function defines what the robot will do.
14      It's called automatically when the protocol runs.
15      """
16
17      # 1. Load labware (define what equipment is on the deck)
18      # Slot numbers 1-11 refer to positions on the robot deck
19      plate = protocol.load_labware('corning_96_wellplate_360ul_flat', '1')
20      tiprack = protocol.load_labware('opentrons_96_tiprack_300ul', '11')
21      reservoir = protocol.load_labware('opentrons_10_tuberack_falcon_4x50ml_6x15ml_conical', '3')
22
23      # 2. Load pipette (define which pipette to use)
24      # 'right' or 'left' refers to the pipette mount position
25      pipette = protocol.load_instrument(
26          'p300_single_gen2',
27          'left',
28          tip_racks=[tiprack]
29      )
30
31      # 3. Protocol steps (define the liquid handling operations)
32
33      # Pick up a tip
34      # pipette.pick_up_tip()
35      pipette.pick_up_tip(tiprack['A1'])
36
37      # Aspirate 20µL from reservoir well A1
38      pipette.aspirate(20, reservoir['A3'])
39
40      # Dispense 20µL to plate well A1
41      pipette.dispense(20, plate['A1'])
42
43      # Repeat for multiple wells
44      for well in ['A2', 'A3', 'A4', 'A5', 'A6']:
45          pipette.aspirate(20, reservoir['A3'])
46          pipette.dispense(20, plate[well])
47
48      # Drop the tip when finished
49      pipette.drop_tip()
```



Everything within run() will be executed


Define and load the labwares

Define and load the modules


Protocol steps

Ref: <https://docs.opentrons.com/python-api/>

# Core Workflow Overview of OT2



# Manual run via Apps




# Communication protocol using Rest API



- The interaction between your computer (the client) and the OT-2 (the server) is managed through a REST API.
- Scripts send HTTP requests to specific URLs on the robot, and the robot responds with JSON-formatted data.

# Automated run via RestAPI--Code walkthrough



The screenshot displays a Jupyter Notebook interface. At the top, the title bar reads 'jupyter HTTP with Image processing NTU\_final\_20260909 Last Checkpoint: 5 hours ago (autosaved)'. Below the title bar is a menu bar with options: File, Edit, View, Insert, Cell, Kernel, Widgets, Help. A status bar indicates 'Not Trusted' and 'Python 3 (ipykernel)'. The main content area shows the tutorial title 'Opentron OT-2 Robot Control and Image Processing Tutorial' and an 'Overview' section. The overview text states: 'This Jupyter notebook provides a comprehensive tutorial for undergraduate students to learn how to use the Opentron OT-2 liquid handling robot. The tutorial covers:'. A numbered list follows: 1. **Basic Opentron Concepts:** Understanding labware, pipettes, and protocol structure; 2. **Robot Connectivity:** How to connect and communicate with the OT-2 robot; 3. **Protocol Development:** Creating and uploading Python protocols; 4. **Real-time Monitoring:** Tracking protocol execution and liquid handling operations; 5. **Image Processing:** Using the built-in camera for automated image capture and analysis; 6. **Data Analysis:** Processing ROI (Region of Interest) data for quantitative analysis. Below the overview is a section titled 'What is the Opentron OT-2?' with the text: 'The Opentron OT-2 is an open-source liquid handling robot designed for laboratory automation. It can:'. A bulleted list follows: • Pipette liquids with high precision (0.5-1000 µL); • Work with standard lab equipment (tips, plates, tubes); • Execute complex protocols automatically; • Capture images during experiments for analysis. The final section is 'Learning Objectives' with the text: 'By the end of this tutorial, you will understand how to:'. A bulleted list follows: • Set up communication with the OT-2 robot; • Define labware and pipette configurations; • Write and upload protocols to the robot; • Monitor protocol execution in real-time; • Process images captured during experiments; • Analyze experimental data using Python.

**Opentron OT-2 Robot Control and Image Processing Tutorial**

**Overview**

This Jupyter notebook provides a comprehensive tutorial for undergraduate students to learn how to use the Opentron OT-2 liquid handling robot. The tutorial covers:

1. **Basic Opentron Concepts:** Understanding labware, pipettes, and protocol structure
2. **Robot Connectivity:** How to connect and communicate with the OT-2 robot
3. **Protocol Development:** Creating and uploading Python protocols
4. **Real-time Monitoring:** Tracking protocol execution and liquid handling operations
5. **Image Processing:** Using the built-in camera for automated image capture and analysis
6. **Data Analysis:** Processing ROI (Region of Interest) data for quantitative analysis

**What is the Opentron OT-2?**

The Opentron OT-2 is an open-source liquid handling robot designed for laboratory automation. It can:

- Pipette liquids with high precision (0.5-1000 µL)
- Work with standard lab equipment (tips, plates, tubes)
- Execute complex protocols automatically
- Capture images during experiments for analysis

**Learning Objectives**




By the end of this tutorial, you will understand how to:

- Set up communication with the OT-2 robot
- Define labware and pipette configurations
- Write and upload protocols to the robot
- Monitor protocol execution in real-time
- Process images captured during experiments
- Analyze experimental data using Python



<https://github.com/MH-BY/MS4672/>

# Automated colour mixing using Opentrons




- Amount of colour x transferred to each well
- Efficiency in liquid transferring



# In-class Project (45mins coding + 15mins run)

(Follow the code provided to achieve automated color mixing)

1. Modify the code provided (protocol\_modified.py) to achieve automated rainbow color mixing in two rows.




2. Name it <matric\_no>.py (e.g. G2201923G.py). Submit your code to:  
<https://github.com/MH-BY/MS4672/tree/main/week-5/Proj-submission>

# **2nd half**


## **OpenCV image manipulation**




# A lot of materials data is visual



Optical microscopy




SEM




Photos of the samples

- Morphology
- Composition
- Homogeneity
- Concentration
- Hardness
- Particle size
- Spectral Absorption




AFM



2D mapping

Image analysis software – Fiji, Gwiddiyon, Ilastik



<https://science-saved.com/afm>  
<https://ieeexplore.ieee.org/abstract/document/10343888>  
<https://chemistry-europe.onlinelibrary.wiley.com/doi/epdf/10.1002/cssc.202402499>



# OpenCV Library



- Image manipulation
  - Read and display
  - Color space and channel manipulation
  - Resize, rescale, cropping
  - Transformations
  - Draw and Masking
- Object and Face detection
  - Contour, thresholding, and edge detection
  - Basic item detection and recognition
- Connecting with AI modules



[Open Computer Vision Library](https://opencv.org/)

*pip install opencv-python*  
*Python >=3.6*


# Image



[https://brandonrohrer.com/convert\\_rgb\\_to\\_grayscale.html](https://brandonrohrer.com/convert_rgb_to_grayscale.html)



# Functionalities




Img



Crop

Resize



Mask & rotate



Transform



CVT



Contour / Edge

[https://brandonrohrer.com/convert\\_rgb\\_to\\_grayscale.html](https://brandonrohrer.com/convert_rgb_to_grayscale.html)



# Code walkthrough

